

MSH-300i digital magnetic stirrer hotplate

2
year warranty

- FUNCTIONAL** › External temperature probe for more accurate temperature control*
- CONVENIENT** › LCD 2 line display of set and actual values
- GOOD VALUE** › Excellent specification for the price, includes stirrer bar and retort stand as standard
- RELIABLE** › 2 year warranty, Grant renowned technical and service support
- VERSATILE** › Large capacity stirring and heating, up to 20L



* Optional

Specifications

	MSH-300i 3.3 kg l: 100 mm w: 190 mm h: 275 mm	 
Mixing speed control range	250 to 1250 rpm	
Maximum stirring volume	20 litres	
Temperature range	+30 to 330°C	
Temperature uniformity	± 3°C	
Plate diameter	160 mm	
Plate material	aluminium alloy	
Retort stand height	320 mm	
Maximum length of stirring bar	50 mm	
Heating power	550 W	

Applications

Clinical/Healthcare - thawing/heating/stirring samples
 Pharmaceutical - heating and stirring samples
 Science education in schools/universities - practical science demonstration and experimentation
 Industrial - QC testing, sample preparation
 Reaction chemistry

Further product information:

Call us:
+ 44 (0) 1763 264 741

Email:
salesdesk@grantinstruments.com