


Freezer racks and boxes

Name	Model	Specification	Dimension (W*D*H)mm	Capacity	Material
Racks	R255A	For 2 inch box, side type;5 layer*5 row	137*692*292	25 units boxes 2"	SUS 304
Racks	R255B	For 2 inch box, drawer type;5 layer*5 row	142*693*290	25 units boxes 2"	
Racks	R335A	For 3 inch box, side type;3 layer*5 row	137*692*260	15 units boxes 3"	
Racks	R335B	For 3 inch box, drawer type;3 layer*5 row	142*693*262	15 units boxes 3"	
Racks	R254A	For 2 inch box, side type;5 layer*4row	137*555*292	20 units boxes 2"	
Racks	R254B	For 2 inch box, drawer type;5 layer*4 row	142*556*290	20 units boxes 2"	
Racks	R244A	For 2 inch box, side type;4 layer*4 row	137*555*234	16 units boxes 2"	
Racks	R244B	For 2 inch box, drawer type;4 layer*4 row	142*556*233	16 units boxes 2"	
Racks	R334A	For 3 inch box, side type;3 layer*4row	137*555*260	12 units boxes 3"	
Racks	R334B	For 3 inch box, drawer type;3 layer*4 row	142*556*262	12 units boxes 3"	
Racks	R253A	For 2 inch box, side type;5 layer*3row	137*418*292	15 units boxes 2"	
Racks	R253B	For 2 inch box, drawer type;5 layer*3 row	142*419*290	15 units boxes 2"	
Racks	R11	For 2 inch box, hike type;11 layer	140*147*630	11 units boxes 2"	
Racks	R8	For 2 inch box, hike type;8 layer	140*147*462	8 units boxes 2"	
Racks	R10	For 2 inch box, hike type;10 layer	140*147*574	10units boxes 2"	



Material of Box:ABS
Size:2"



Material of Box:Paper Characteristics
Size:3"

Note:

"A" means side type rack; "B" means drawer type rack;